




## Topi-Click Use Instructions & Picture Guide




**STEP 1:**  
Lift off the cap of the Topi-CLICK® container by placing your thumb under the tab on the protective outer cap and pressing upward.



**STEP 2:**  
Turn the base to dispense the number of CLICKs specified in the instructions on your prescription label.




**STEP 3:**  
Watch the visual indicators when rotating base. Be sure to leave the base in the registered CLICK position by aligning the base column with the tube node after each use.




**STEP 4:**  
Apply topical medication by rubbing the area of skin as directed by your pharmacist using the built-in applicator. Continue rubbing until topical is completely absorbed into the skin.

**IMPORTANT:** The built-in applicator lessens the risk of cross-contamination and ensures the full dose is applied to the intended area.



**STEP 5:**  
Replace cap after each dose to prevent evaporation. Be sure cap snaps on tightly. **Note** Be sure to align ports when using the Topi-CLICK® 35 3Port™ cap and applicator.



**STEP 6:**  
The bottom of the plunger is the indicating line; when it reaches the embossed line below the word "Refill" approximately 32 CLICKs are left. Approximately 12 CLICKs remain at the line above the word "Refill".